
Quick Start Guide



Immuneservice's WST-1 Cell Proliferation Kit (ready-to-use) provides an easy, fast and accurate method to determine cell proliferation, cell viability and cytotoxicity.

5-Step Protocol:

1. Culture $0.1-5 \times 10^4$ cells/well in 96-well microtiter plate in a final volume of 100 μ l/well culture medium for 24-96 hours.
2. Add 10 μ l/well of WST-1 ready-to-use Cell Proliferation Kit.
3. Incubate the microtiter plate for 1-4 hours.
4. Measure absorbance at 420 - 480 nm, with a reference wavelength above 600 nm.
5. Analyze data.

If you have any questions, please refer to our detailed user manual or contact us by email: support@immuneservice.com.

For research use only. Not intended for diagnostic or therapeutic procedures.
