

---

## Quick Start Guide

---



Immunservice's WST-1 CTLL-2 Cell Proliferation Kit (ready-to-use) provides an easy, fast and accurate method to determine cell proliferation, cell viability and cytotoxicity with the cytotoxic murine T-cell line CTLL-2.

### 5-Step Protocol:

1. Culture  $4 \times 10^3$  CTLL-2 cells/well on 96-well microtiter plate in a final volume of 100  $\mu$ l/well culture medium.
2. Add 10  $\mu$ l/well of WST-1 CTLL-2 Cell Proliferation Kit.
3. Incubate the microtiter plate for 1–4 hours.
4. Measure absorbance at 420 – 480 nm, with a reference wavelength above 600 nm.
5. Analyze data.

If you have any questions, please refer to our detailed user manual or contact us by email: [support@immunservice.com](mailto:support@immunservice.com).

---

For research use only. Not intended for diagnostic or therapeutic procedures.

---